

Product Information


DSPE PEG reagent, DSPE-PEG-folate, MW 1,000, Purity $\geq 95\%$

Cat. No.: X25-10-YM197

Size: 25 mg; 50 mg; 100 mg; 200 mg

Synonym: DSPE-PEG-folate

This product is for research use only and is not intended for diagnostic use.



Product Information

Description	Stearic acid-PEG-FITC comprises the functional group FITC, with reactivity marked by amine, thiol, and alcohol.
Functional Group 1	FITC
Reactive Group 1	Amine
Reactive Group 2	Thiol
Reactive Group 3	Alcohol
Form	Solid
Purity	$\geq 95\%$
Impurities	No visible impurities to the naked eye.
Identity	Confirmed by NMR.
Applications	Stearic acid-PEG-FITC plays a key role in studying membrane interactions and dynamics due to the presence of stearic acid, which facilitates membrane incorporation.
Storage	Store at -20°C .